

### Anti-Bcl2-L-13 antibody ab27795

1 References 画像数 1

医薬用外劇物

#### 製品の概要

製品名	Anti-Bcl2-L-13 antibody
製品の詳細	Rabbit polyclonal to Bcl2-L-13
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to Mouse Bcl2-L-13. (Peptide available as <a href="#">ab41715</a> )
ポジティブ・コントロール	Jurkat cell lysate.
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3671 Bcl-Rambo Antibody. 3671-100 is the same size as the 100 µg size of ab27795.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Thimerosal (merthiolate) Constituents: PBS, 20% Glycerol (glycerin, glycerine), 0.1% BSA
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab27795の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 4 µg/ml. Predicted molecular weight: 53 kDa.

## ターゲット情報

機能	May promote the activation of caspase-3 and apoptosis.
組織特異性	Ubiquitous, with the highest levels of expression in heart, placenta and pancreas.
配列類似性	Belongs to the Bcl-2 family.
細胞内局在	Nucleus and Mitochondrion membrane. Nucleus.

## 画像



Anti-Bcl2-L-13 antibody (ab27795) at 4 µg/ml + Jurkat cell lysate

**Predicted band size:** 53 kDa

**Observed band size:** 56 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors