

Biotin Anti-beta-casomorphin antibody ab47859

医薬用外劇物

製品の概要

製品名	Biotin Anti-beta-casomorphin antibody
製品の詳細	Biotin Rabbit polyclonal to beta-casomorphin
由来種	Rabbit
標識	Biotin
アプリケーション	適用あり: ELISA, RIA
種交差性	交差種: Cow
免疫原	Synthetic peptide corresponding to Cow beta-casomorphin conjugated to keyhole limpet haemocyanin. Sequence: CYPFPGPI

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.50 Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 50% Glycerol
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab47859の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.

ターゲット情報

関連性 β -casomorphin (β -CM) is an exogenous bioactive peptide derived from casein, a major protein in milk and milk products, which has opioid activity. The circulation of this peptide into children's immature central nervous system is believed to be related with the "Sudden Infant Death Syndrome (SIDS)". It is believed that this peptide produces an inhibitory effect in the respiratory centre in the brain stem, and that this leads to apnea and death.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors