

### Biotin Anti-Borrelia burgdorferi antibody ab69158

#### 製品の概要

---

製品名	Biotin Anti-Borrelia burgdorferi antibody
製品の詳細	Biotin Rabbit polyclonal to Borrelia burgdorferi
由来種	Rabbit
標識	Biotin
特異性	ab69158 reacts with Borrelia burgdorferi in bacterial and infected tissue samples.
アプリケーション	<b>適用あり:</b> IHC, ICC/IF
種交差性	<b>交差種:</b> Borrelia burgdorferi
免疫原	Tissue, cells or virus corresponding to Borrelia burgdorferi. Borrelia burgdorferi whole cell preparation.

#### 特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

---

The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab69158の使用に適用されます**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent dilution.

## ターゲット情報

---

**関連性** *Borrelia burgdorferi* is a spirochete and the cause of Lyme disease, a tick transmitted illness of humans and animals. *B. burgdorferi* may persist in humans and animals for months or years following initial infection, despite a robust humoral immune response. *B. burgdorferi* resembles other spirochetes in that it is a highly specialized, motile, two-membrane, spiral shaped bacteria which lives primarily as an extracellular pathogen. *B. burgdorferi* has an unusual genome compared with other eubacteria which includes a linear chromosome approximately one megabase in size and numerous linear and circular plasmids.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors