

Biotin Anti-CD200 / OX2 antibody [MRC OX90] ab33735

製品の概要

製品名	Biotin Anti-CD200 / OX2 antibody [MRC OX90]
製品の詳細	Biotin Rat monoclonal [MRC OX90] to CD200 / OX2
由来種	Rat
標識	Biotin
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Fusion protein corresponding to Mouse CD200/ OX2. Database link: O54901

特記事項

Spleen cells from immunised rats were fused with cells of the rat Y3 myeloma cell line.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
精製度	Protein G purified
特記事項(精製)	This purified antibody was prepared by affinity chromatography on Protein G from tissue culture supernatant.
一次抗体 備考	Spleen cells from immunised rats were fused with cells of the rat Y3 myeloma cell line.
ポリモノ	モノクローナル
クローン名	MRC OX90

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab33735の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl.

ターゲット情報

機能	Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.
配列類似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors