

Biotin Anti-E. coli antibody ab20640

24 References

製品の概要

製品名	Biotin Anti-E. coli antibody
製品の詳細	Biotin Rabbit polyclonal to E. coli
由来種	Rabbit
標識	Biotin
特異性	Many "O" and "K" antigenic serotypes of Escherichia coli. Will remove E. coli proteins from recombinant preparations. Tested specifically with the strains O18, O20, O44, O55, O111, O112, O125, O157, K12.
アプリケーション	適用あり: ELISA, ICC/IF
種交差性	交差種: Escherichia coli
免疫原	Tissue, cells or virus corresponding to Escherichia coli E. coli. Mixture of E. coli serotypes: K12, O111, O55, O125, O20, and O157
特記事項	<p>Covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Protein A purified
特記事項 (精製)	Protein A chromatography.

ポリ/モノ
アイソタイプ

ポリクローナル
IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab20640の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.

ターゲット情報

関連性 Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriae - bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, especially where there is an immature immune system or immunosuppression, or where the subtype of organism has acquired the ability to produce pathogenic toxins.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors