

### Biotin Anti-GM-CSF antibody ab84259

#### 製品の概要

製品名	Biotin Anti-GM-CSF antibody
製品の詳細	Biotin Rabbit polyclonal to GM-CSF
由来種	Rabbit
標識	Biotin
アプリケーション	<b>適用あり:</b> Sandwich ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Highly pure (>98%) recombinant human GM-CSF.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Lyophilized:Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituent: PBS
精製度	Protein A purified
特記事項(精製)	ab84259 is sterile filtered.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab84259の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use at an assay dependent concentration. Allows the detection of at least 0.2 - 0.4 ng/well of recombinant human GM-CSF. Can be paired for Sandwich ELISA with <b>ab9667</b> . Can be used as detection antibody when paired with recommended pair.
WB		Use at an assay dependent concentration. Predicted molecular weight: 16 kDa. The detection limit for recombinant human GM-CSF is 1.5 - 3.0 ng/lane under either reducing or non-reducing conditions.

#### ターゲット情報

機能	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.
配列類似性	Belongs to the GM-CSF family.
細胞内局在	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors