

Biotin Anti-ICAM1 antibody ab7815

★★★★★ [3 Abreviews](#) [7 References](#) [画像数 1](#)

製品の概要

製品名	Biotin Anti-ICAM1 antibody
製品の詳細	Biotin Rabbit polyclonal to ICAM1
由来種	Rabbit
標識	Biotin
アプリケーション	適用あり: Sandwich ELISA
種交差性	交差種: Human 交差が予測される動物種: Chimpanzee, Baboon, Cynomolgus monkey, Rhesus monkey 
免疫原	Full length protein corresponding to Human ICAM1. Database link: P05362
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.065% Sodium azide Constituent: 0.1% BSA
特記事項(精製)	0.2 µm filtered
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab7815の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use a concentration of 0.5 µg/ml. Can be paired for Sandwich ELISA with Mouse monoclonal [MEM-111] to ICAM1 (ab2213) . For sandwich ELISA, use this antibody as Detection at 0.5µg/ml with ab2213 as Capture.

ターゲット情報

機能	ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). During leukocyte trans-endothelial migration, ICAM1 engagement promotes the assembly of endothelial apical cups through ARHGEF26/SGEF and RHO G activation. In case of rhinovirus infection acts as a cellular receptor for the virus.
配列類似性	Belongs to the immunoglobulin superfamily. ICAM family. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
翻訳後修飾	Monoubiquitinated, which is promoted by MARCH9 and leads to endocytosis.
細胞内局在	Membrane.

画像



Standard Curve for ICAM1 (Analyte: **ab82125**) dilution range 1pg/ml to 1ug/ml using Capture Antibody Mouse monoclonal [MEM-111] to ICAM1 (**ab2213**) at 5ug/ml and Detector Antibody Rabbit polyclonal to ICAM1 (Biotin) (ab7815) at 0.5ug/ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors