

### Biotin Anti-Melanoma Inhibitory Activity antibody ab84235

★★★★★ **1 Abreviews**

#### 製品の概要

製品名	Biotin Anti-Melanoma Inhibitory Activity antibody
製品の詳細	Biotin Rabbit polyclonal to Melanoma Inhibitory Activity
由来種	Rabbit
標識	Biotin
アプリケーション	<b>適用あり:</b> Sandwich ELISA, Flow Cyt, WB
種交差性	<b>交差種:</b> Human
免疫原	Highly pure (>98%) recombinant full length Human Melanoma Inhibitory Activity
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Lyophilized:Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituent: PBS
特記事項(精製)	ab84235 was purified by affinity chromatography and then biotinylated.Filter sterilised.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、**次のテスト済みアプリケーションにおけるab84235の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use at an assay dependent concentration.
Flow Cyt	★★★★★ (1)	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

## ターゲット情報

機能	Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.
組織特異性	All malignant melanoma cell lines tested and infrequently in glioma cell lines.
配列類似性	Belongs to the MIA/OTOR family. Contains 1 SH3 domain.
翻訳後修飾	May possess two intramolecular disulfide bonds.
細胞内局在	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors