

### Biotin Anti-PDGF B antibody ab84265

#### 製品の概要

製品名	Biotin Anti-PDGF B antibody
製品の詳細	Biotin Rabbit polyclonal to PDGF B
由来種	Rabbit
標識	Biotin
アプリケーション	<b>適用あり:</b> Sandwich ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Highly pure (>98%) recombinant human PDGF BB
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Lyophilized:Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituent: PBS
特記事項(精製)	ab84265 is affinity purified and sterile filtered.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab84265の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
<b>Sandwich ELISA</b>		Use at an assay dependent concentration. To be used with <b>ab9704</b> as capture antibody and <b>ab9706</b> as protein standard. Allows the detection of at least 0.2 - 0.4 ng/well of recombinant human PDGF BB.
<b>WB</b>		Use at an assay dependent concentration. Predicted molecular weight: 27 kDa. The detection limit for recombinant human PDGF BB is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.

## ターゲット情報

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors