

# Anti-Blood Group H ab antigen antibody [97-I] ab24213

★★★★☆ [1 Abreviews](#) [1 References](#)

### 製品の概要

---

製品名	Anti-Blood Group H ab antigen antibody [97-I]
製品の詳細	Mouse monoclonal [97-I] to Blood Group H ab antigen
由来種	Mouse
特異性	H antigen of red blood cells and body fluids.
アプリケーション	<b>適用あり:</b> ELISA, IP, Flow Cyt, WB
種交差性	<b>交差種:</b> Human
免疫原	Tissue, cells or virus corresponding to Human Blood Group H ab antigen. H antigen of red blood cells and body fluids.
エピトープ	A carbohydrate moiety (the exact composition of the epitope is not mapped).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	97-I
アイソタイプ	IgM

### アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab24213の使用に適用されず**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. <b>ab91545</b> - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody.
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors