

Anti-Bordetella pertussis toxin/ptxD antibody [56.2] ab37686

1 References

製品の概要

製品名	Anti-Bordetella pertussis toxin/ptxD antibody [56.2]
製品の詳細	Mouse monoclonal [56.2] to Bordetella pertussis toxin/ptxD
由来種	Mouse
特異性	This product reacts with different epitope than ab37574 , ab37684 and ab37687 .
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Bordetella pertussis
免疫原	<i>Bordetela pertussis</i> Toxin B-Oligomer.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.40 Constituents: 2.9% Sodium chloride, 0.0268% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	56.2
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab37686の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Detects a band of approximately 11.7 kDa.

ターゲット情報

関連性 Pertussis toxin contains five different chains, S1-S5. They are organized into 2 functional subunits: A, composed of S1 (which is toxic) and B, containing S2, S3, S5, and two copies of S4 (B binds to the membrane receptors). Dimers of S2-S4 and S3-S4 are held together by S5

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors