

Anti-BOTOX antibody ab20786

★★★★☆ 1 Abreviews

製品の概要

製品名	Anti-BOTOX antibody
製品の詳細	Rabbit polyclonal to BOTOX
由来種	Rabbit
特異性	This antibody is specific for Clostridium botulinum Type A Neurotoxin Light Chain.
アプリケーション	適用あり: Dot blot, ELISA
種交差性	交差種: Species independent
免疫原	Other Immunogen Type corresponding to BOTOX. Type A Neurotoxin Light Chain
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 2.68% PBS, 20% Glycerol (glycerin, glycerine)
精製度	Protein G purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab20786の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

ターゲット情報

関連性 Clostridium botulinum neurotoxins are zinc-dependent proteases that block the release of acetylcholine from peripheral cholinergic nerve endings by targeting specific membrane proteins including synaptosomal-associated proteins (SNAP), vesicle-associated membrane proteins (VAMP), and syntaxins. Neurotoxin A is characteristically associated with human botulism and known to target and hydrolyse the 197-Gln-|-Arg-198 bond of SNAP-25.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors