

# Anti-BRAG2 antibody ab122442

**1 References** [画像数 2](#)

### 製品の概要

製品名	Anti-BRAG2 antibody
製品の詳細	Rabbit polyclonal to BRAG2
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, ICC/IF
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment corresponding to Human BRAG2 aa 28-96 (N terminal). Sequence: PSSETGTSLDSPSAYPQGPLVPGSSLSPDHYEHTSVGAYGLY SGPPGQQQ RTRRPKLQHSTSILRKQAE  Database link: <a href="#">Q6DN90</a>  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
ポジティブ・コントロール	Human skin tissue; Human cell line A431
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee** **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab122442の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.

## ターゲット情報


機能	In addition to accelerate GTP gamma S binding by ARFs of all three classes, it appears to function preferentially as a guanine nucleotide exchange protein for ARF6, mediating internalisation of beta-1 integrin.
組織特異性	Expressed in brain, ovary, heart, lung, liver, kidney and leukocytes. Moderate expression was also detected in lung, skeletal muscle, placenta, small intestine, pancreas, spleen and testis.
配列類似性	Belongs to the BRAG family. Contains 1 IQ domain. Contains 1 PH domain. Contains 1 SEC7 domain.
細胞内局在	Cytoplasm. Nucleus. At steady state, may be preferentially cytosolic.

## 画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-BRAG2 antibody (ab122442)

ab122442, at a 1/50 dilution, staining BRAG2 in paraffin embedded Human skin tissue by Immunohistochemistry.



Immunocytochemistry/ Immunofluorescence - Anti-BRAG2 antibody (ab122442)

ab122442, at 4  $\mu\text{g/ml}$ , staining BRAG2 in Human A431 cell line, in green, by Immunofluorescence.

Fixation/Permeabilization: PFA/Triton X-100.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors