

Anti-BTG2/PC3 antibody ab85051

★★★★★ 1 Abreviews 8 References 画像数 2

医薬用外毒物

製品の概要

製品名	Anti-BTG2/PC3 antibody
製品の詳細	Rabbit polyclonal to BTG2/PC3
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human BTG2/PC3 (C terminal). A 15 amino acid synthetic peptide from the C-terminal residues. Database link: NM_006763
ポジティブ・コントロール	WB: HepG2 whole cell lysate (ab7900). IHC-P: Human fetal colon tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute in 200ul sterile H2O.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: 2% BSA
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab85051の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB	★★★★☆ (1)	1/200 - 1/1000. Detects a band of approximately 18 kDa (predicted molecular weight: 18 kDa).
IHC-P		1/100 - 1/500.

ターゲット情報


機能	Involved in cell cycle regulation. Could be involved in the growth arrest and differentiation of the neuronal precursors (By similarity). Anti-proliferative protein. Modulates transcription regulation mediated by ESR1. Involved in neurite outgrowth.
配列類似性	Belongs to the BTG family.

画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human fetal colon tissue labelling BTG2/PC3 with ab85051 at a dilution of 1/100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-BTG2/PC3 antibody (ab85051)



Anti-BTG2/PC3 antibody (ab85051) at 1/250 dilution + HepG2 cell lysates

Predicted band size: 18 kDa

Observed band size: 18 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors