

Anti-C Reactive Protein antibody [EPR283Y] ab51016

リコンビナント **RabMAb**

画像数 2

製品の概要

製品名	Anti-C Reactive Protein antibody [EPR283Y]
製品の詳細	Rabbit monoclonal [EPR283Y] to C Reactive Protein
由来種	Rabbit
アプリケーション	適用あり: WB 適用なし: Flow Cyt or IP
種交差性	交差種: Human
免疫原	Synthetic peptide within Human C Reactive Protein aa 1-100 (N terminal). The exact sequence is proprietary.
ポジティブ・コントロール	Human serum
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル

クローン名 EPR283Y
アイソタイプ IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab51016の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000. Detects a band of approximately 25 kDa (predicted molecular weight: 25 kDa).

追加情報 Is unsuitable for Flow Cyt or IP.

ターゲット情報


機能 Displays several functions associated with host defense: it promotes agglutination, bacterial capsular swelling, phagocytosis and complement fixation through its calcium-dependent binding to phosphorylcholine. Can interact with DNA and histones and may scavenge nuclear material released from damaged circulating cells.

組織特異性 Found in plasma.

配列類似性 Belongs to the pentaxin family.
Contains 1 pentaxin domain.

細胞内局在 Secreted.

画像



Anti-C Reactive Protein antibody [EPR283Y] (ab51016) at 1/1000 dilution + Human serum at 10 µg

Secondary

Goat anti-rabbit HRP labeled at 1/2000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-C Reactive Protein antibody [EPR283Y]
(ab51016)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors