

Anti-CALML5 antibody ab122665

画像数 2

製品の概要

製品名	Anti-CALML5 antibody
製品の詳細	Rabbit polyclonal to CALML5
由来種	Rabbit
アプリケーション	適用あり: IHC-P, ICC/IF
種交差性	交差種: Human
免疫原	antigen sequence: GELTPEEEAQ YKKAFAVDV DGNGTINAQE LGAALKATGK NLSEAQLRKL ISEVDGDDG EISFQEFLTA, corresponding to N terminal amino acids 3-72 of Human CALML5.

[Run BLAST with Expasy](#) [Run BLAST with NCBI](#)

ポジティブ・コントロール Human skin tissue. CALML5 over-expression lysate.

特記事項 The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab122665の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/2500 - 1/5000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.

ターゲット情報


機能	Binds calcium. May be involved in terminal differentiation of keratinocytes.
組織特異性	Particularly abundant in the epidermis where its expression is directly related to keratinocyte differentiation. Very low expression in lung.
配列類似性	Contains 4 EF-hand domains.

画像



ab122665 at 1/2500 dilution staining CALML5 in paraffin-embedded Human skin tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CALML5 antibody (ab122665)



Immunocytochemistry/ Immunofluorescence - Anti-CALML5 antibody (ab122665)

Immunofluorescent staining of Human cell line A-431 shows positivity in plasma membrane and cytoplasm. Recommended concentration of ab122665 1-4 $\mu\text{g/ml}$. Cells treated with PFA/Triton X-100.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors