


Anti-Calponin 1 antibody [EP798Y] ab46794

リコンビナント **RabMAb**

★★★★★ **20 Abreviews** **234 References** 画像数 21


製品の概要

製品名	Anti-Calponin 1 antibody [EP798Y]
製品の詳細	Rabbit monoclonal [EP798Y] to Calponin 1
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P, ICC/IF 適用なし: Flow Cyt
種交差性	交差種: Mouse, Rat, Human, Pig 交差が予測される動物種: Sheep 
免疫原	Synthetic peptide within Human Calponin 1 aa 250-350 (C terminal). The exact sequence is proprietary. Database link: P51911
ポジティブ・コントロール	ICC/IF: C2C12 and HeLa cells; Pig and mouse aortic smooth muscle cells. IHC-P: Rat lung tissue; Mouse cardiac muscle tissue; Human lung carcinoma, kidney, lung, tonsil, uterus, smooth muscle and skeletal muscle tissues. WB: Human bladder lysates, Pig heart lysates, Mouse bladder lysates and Rat bladder lysates.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
解離定数 (K _D 値)	K _D = 1.73 x 10 ⁻¹⁰ M

10⁻¹⁰



[Learn more about K_D](#)

バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol, 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EP798Y
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab46794の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB	★★★★★ (9)	1/5000. Detects a band of approximately 34 kDa. For unpurified use at 1/20,000.
IHC-P	★★★★★ (3)	1/1000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. For unpurified use at 1/100 - 1/250.
ICC/IF	★★★★★ (4)	1/500. For unpurified use at 1/100 - 1/250.

追加情報 Is unsuitable for Flow Cyt.

ターゲット情報


機能	Thin filament-associated protein that is implicated in the regulation and modulation of smooth muscle contraction. It is capable of binding to actin, calmodulin, troponin C and tropomyosin. The interaction of calponin with actin inhibits the actomyosin Mg-ATPase activity.
組織特異性	Smooth muscle, and tissues containing significant amounts of smooth muscle.
配列類似性	Belongs to the calponin family. Contains 3 calponin-like repeats. Contains 1 CH (calponin-homology) domain.

画像



Immunocytochemistry/ Immunofluorescence - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Immunocytochemistry/ Immunofluorescence analysis of C2C12 (Mouse myoblasts myoblast) cells labeling Calponin 1 with purified ab46794 at 1/500 dilution. Cells were fixed in 4% paraformaldehyde and permeabilized with 0.1% Triton X-100. Cells were counterstained with **ab195889** Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) 1/200 (2.5 µg/ml). **ab150077** Goat anti rabbit IgG(Alexa Fluor® 488) was used as the secondary antibody at 1:1000 dilution. DAPI nuclear counterstain. PBS instead of the primary antibody was used as the secondary antibody only control.



Western blot - Anti-Calponin 1 antibody [EP798Y] (ab46794)

All lanes : Anti-Calponin 1 antibody [EP798Y] (ab46794) at 1/50000 dilution (purified)

Lane 1 : Human bladder lysates

Lane 2 : Mouse bladder lysates


Lane 3 : Rat bladder lysates

Lysates/proteins at 15 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG (HRP) with minimal cross-reactivity with human IgG at 1/2000 dilution


Blocking and diluting buffer: 5% NFDM/TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human lung carcinoma tissue sections labeling Calponin 1 with purified ab46794 at a 1:1000 dilution. Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0). Tissue was counterstained with hematoxylin. ImmunoHistoProbe one step HRP Polymer (ready to use) secondary antibody was used.

PBS instead of the primary antibody was used as the negative control (inset).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)



All lanes : Anti-Calponin 1 antibody [EP798Y] (ab46794)


Lanes 1-2 : Non-Failing pig LV-myocardium tissue lysate.

Lanes 3-4 : Failing pig LV-myocardium tissue lysate injected with DOX.

Ab46794 was used to stain Calponin 1 in Neonatal piglets injected with Dox (Failing LV myocardium) and PBS (Non-Failing LV myocardium). At the protein level a more than twofold increase in Calponin 1 was observed in Dox-injected animals compared to controls (PBS).

Western blot - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Torrado Met al. Targeted gene-silencing reveals the functional significance of myocardin signaling in the failing heart. PLoS One 6:e26392 (2011). Reproduced under the Creative Commons license <http://creativecommons.org/licenses/by/4.0/>




Immunocytochemistry/ Immunofluorescence - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Paraformaldehyde-fixed, 0.25% Triton X-100 permeabilized mouse thoracic aortic smooth muscle cells labeling Calponin 1 using ab46794 at 1/100 dilution in ICC/IF, followed by a Goat Anti-Rabbit IgG H&L (Alexa Fluor 488) ([ab150077](#)) at 1/400 dilution.

1.5% BSA used as blocking agent for 30 minutes at 25°C. Incubated with primary antibody for 24 hours at 4°C.

VSMCs were seeded to 35-mm plates in a low density avoiding overlapping of cells. After fixation, VSMCs were treated with 0.25% Triton X-100 for 20 minutes.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of mouse cardiac muscle tissue sections labeling Calponin 1 with purified ab46794 at 1:1000 dilution. Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0). Tissue was counterstained with hematoxylin.

ImmunoHistoProbe one step HRP Polymer (ready to use) secondary antibody was used.

PBS instead of the primary antibody was used as the negative control (inset).




Western blot - Anti-Calponin 1 antibody [EP798Y]
(ab46794)

Anti-Calponin 1 antibody [EP798Y] (ab46794) at 1/5000 dilution
(purified) + Pig heart lysates at 15 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


Blocking and diluting buffer: 5% NFDM/TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-Calponin 1 antibody
[EP798Y] (ab46794)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of rat lung tissue sections labeling Calponin 1 with purified ab46794 at 1:1000 dilution. Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0). Tissue was counterstained with hematoxylin. ImmunoHistoProbe one step HRP Polymer (ready to use) secondary antibody.

PBS instead of the primary antibody was used as the negative control (inset).




Immunocytochemistry/ Immunofluorescence - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Image kindly supplied by Mrs Jordan Carbary through Abreview

Unpurified ab46794 staining Calponin 1 in porcine aortic smooth muscle cells by Immunocytochemistry/ Immunofluorescence.

The cells were paraformaldehyde fixed, permeabilized in 0.1% Triton X-100. Samples were then incubated with primary antibody at 1/50 for 1 hour at 25°C. The secondary antibody used was **ab6717** Goat polyclonal to Rabbit IgG - H&L (FITC) (green) used at a 1/400 dilution.




Immunocytochemistry/ Immunofluorescence - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Herington et al PLoS One. 2015 Nov 24;10(11):e0143243. doi: 10.1371/journal.pone.0143243. eCollection 2015. Fig 1. Reproduced under the Creative Commons license <http://creativecommons.org/licenses/by/4.0/>

Representative photomicrograph of mouse UT-myo cells (Left panel) and murine uterine myometrium (Right panel) stained with smooth muscle cell markers, alpha-SMA (red) and ab46794 (green) and DAPI (blue). UT-myo cells and whole-mount uterine tissue were collected from day 19 of mouse pregnancy. The placenta and embryo were removed from whole-mount tissue sections.

For full details please see paper.




OIR-D Scanning - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Equilibrium disassociation constant (K_D)

Learn more about K_D

[Click here to learn more about \$K_D\$](#)




Western blot - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Anti-Calponin 1 antibody [EP798Y] (ab46794) at 1/20000 dilution (unpurified) + Human bladder lysate at 10 µg

Secondary

Goat anti-rabbit HRP at 1/2000 dilution


Observed band size: 34 kDa



Immunocytochemistry/ Immunofluorescence - Anti-Calponin 1 antibody [EP798Y] (ab46794)


ICC/IF image of unpurified ab46794 stained HeLa (Human epithelial cell line from cervix adenocarcinoma) cells.

Cells were 100% methanol fixed (5 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilize the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab46794, 5µg/ml) overnight at +4°C. The secondary antibody (green) was DyLight® 488 goat anti-rabbit IgG - H&L, pre-adsorbed (**ab96899**) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Immunohistochemical staining of paraffin-embedded human smooth muscle using unpurified ab46794 at 1/100 dilution



Western blot - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Image courtesy of an anonymous Abreview.

All lanes : Anti-Calponin 1 antibody [EP798Y] (ab46794) at 1/20000 dilution (unpurified)

All lanes : Whole tissue lysate prepared from bovine lymph vessels

Lysates/proteins at 10 µg per lane.


Secondary

All lanes : HRP conjugated goat anti-rabbit polyclonal at 1/5000 dilution

Developed using the ECL technique.


Observed band size: 34 kDa

Exposure time: 1 minute




Unpurified ab46794 showing positive staining in normal lung vessel tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)




Unpurified ab46794 showing positive staining in normal kidney vessels tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)




Unpurified ab46794 showing positive staining in normal tonsil vessel tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)



Unpurified ab46794 showing positive staining in normal uterus tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)



Unpurified ab46794 showing negative staining in skeletal muscle tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calponin 1 antibody [EP798Y] (ab46794)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Calponin 1 antibody [EP798Y] (ab46794)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors