

# Anti-Calretinin antibody ab702

★★★★☆ **4 Abreviews**   **59 References**   画像数 1

### 製品の概要

製品名	Anti-Calretinin antibody
製品の詳細	Rabbit polyclonal to Calretinin
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Chicken, Cat, Dog, Xenopus laevis 
免疫原	Full length protein corresponding to Calretinin.
特記事項	<p>Contains a proprietary preservative that is not sodium azide or thimerosal.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituents: PBS, Carrier protein  Van Gogh Yellow Diluent
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
ミエローマ	unknown
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab702の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	Use at an assay dependent concentration. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. For an ABC system, dilute 1/50-1/100 and incubate for 30-60 minutes at RT.

## ターゲット情報

機能	Calretinin is a calcium-binding protein which is abundant in auditory neurons.
組織特異性	Brain.
配列類似性	Belongs to the calbindin family. Contains 6 EF-hand domains.

## 画像



Paraffin embedded mesothelioma (human) stained for Calretinin with ab702 (1/100 dilution) in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Calretinin antibody (ab702)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors