

Anti-Caspase-8 antibody ab25901

★★★★☆ [3 Abreviews](#) [136 References](#) [画像数 6](#)

製品の概要

製品名	Anti-Caspase-8 antibody
製品の詳細	Rabbit polyclonal to Caspase-8
由来種	Rabbit
アプリケーション	適用あり: ICC, IHC-P, ICC/IF, WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide corresponding to Human Caspase-8 (C terminal). Synthetic peptide (Human) - a 16 amino acid peptide from near the carboxy-terminus of human Caspase-8 isoform A. Database link: Q14790

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
特記事項(精製)	Caspase-8 Antibody is affinity chromatography purified via peptide column.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee


Abpromise保証は、次のテスト済みアプリケーションにおけるab25901の使用に適用されず

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration. Human samples only
ICC/IF		Use a concentration of 20 µg/ml. Human samples only
WB	★★★★☆ (2)	Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 55 kDa (predicted molecular weight: 55 kDa).


ターゲット情報

機能	<p>Most upstream protease of the activation cascade of caspases responsible for the TNFRSF6/FAS mediated and TNFRSF1A induced cell death. Binding to the adapter molecule FADD recruits it to either receptor. The resulting aggregate called death-inducing signaling complex (DISC) performs CASP8 proteolytic activation. The active dimeric enzyme is then liberated from the DISC and free to activate downstream apoptotic proteases. Proteolytic fragments of the N-terminal propeptide (termed CAP3, CAP5 and CAP6) are likely retained in the DISC. Cleaves and activates CASP3, CASP4, CASP6, CASP7, CASP9 and CASP10. May participate in the GZMB apoptotic pathways. Cleaves ADPRT. Hydrolyzes the small-molecule substrate, Ac-Asp-Glu-Val-Asp-AMC. Likely target for the cowpox virus CRMA death inhibitory protein. Isoform 5, isoform 6, isoform 7 and isoform 8 lack the catalytic site and may interfere with the pro-apoptotic activity of the complex.</p>
組織特異性	<p>Isoform 1, isoform 5 and isoform 7 are expressed in a wide variety of tissues. Highest expression in peripheral blood leukocytes, spleen, thymus and liver. Barely detectable in brain, testis and skeletal muscle.</p>
関連疾患	<p>Defects in CASP8 are the cause of caspase-8 deficiency (CASP8D) [MIM:607271]. CASP8D is a disorder resembling autoimmune lymphoproliferative syndrome (ALPS). It is characterized by lymphadenopathy, splenomegaly, and defective CD95-induced apoptosis of peripheral blood lymphocytes (PBLs). It leads to defects in activation of T-lymphocytes, B-lymphocytes, and natural killer cells leading to immunodeficiency characterized by recurrent sinopulmonary and herpes simplex virus infections and poor responses to immunization.</p>
配列類似性	<p>Belongs to the peptidase C14A family. Contains 2 DED (death effector) domains.</p>
ドメイン	<p>Isoform 9 contains a N-terminal extension that is required for interaction with the BCAP31 complex.</p>
翻訳後修飾	<p>Generation of the subunits requires association with the death-inducing signaling complex (DISC), whereas additional processing is likely due to the autocatalytic activity of the activated protease. GZMB and CASP10 can be involved in these processing events. Phosphorylated upon DNA damage, probably by ATM or ATR.</p>
細胞内局在	<p>Cytoplasm.</p>




Immunocytochemistry/Immunofluorescence of Jurkat cells labeling Caspase-8 with Anti-Caspase-8 antibody (ab25901) at 20µg/ml.

Immunocytochemistry/ Immunofluorescence - Anti-Caspase-8 antibody (ab25901)



Immunohistochemistry of Caspase-8 in human spleen tissue with Caspase-8 antibody at 5 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Caspase-8 antibody (ab25901)



Western blot - Anti-Caspase-8 antibody (ab25901)


All lanes : Anti-Caspase-8 antibody (ab25901) at 0.5 µg/ml

Lane 1 : Human spleen

Lane 2 : Mouse spleen

Lane 3 : Rat spleen

Predicted band size: 55 kDa



Western blot - Anti-Caspase-8 antibody (ab25901)

All lanes : Anti-Caspase-8 antibody (ab25901) at 1 µg/ml

Lane 1 : Jurkat cell lysates with Absence of blocking buffer

Lane 2 : Jurkat cell lysates with Presence of blocking buffer

Predicted band size: 55 kDa



Immunocytochemistry/ Immunofluorescence - Anti-Caspase-8 antibody (ab25901)

ab25901 at 2ug/ml staining Caspase-8 in Jurkat cells by ICC/IF



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Caspase-8 antibody (ab25901)

IHC image of Caspase-8 staining in human spleen formalin fixed paraffin embedded tissue section, performed on a Leica Bond system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with **ab25901**, 1µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors