

Anti-CD11d antibody [CA11.8H2] ab34096

1 References

製品の概要

製品名	Anti-CD11d antibody [CA11.8H2]
製品の詳細	Mouse monoclonal [CA11.8H2] to CD11d
由来種	Mouse
特異性	ab34096 recognises the canine CD11d cell surface antigen (the integrin α sub-unit).
アプリケーション	適用あり: Flow Cyt 適用なし: IHC-P
種交差性	交差種: Dog
免疫原	Full length native protein (purified) corresponding to Dog CD11d. Affinity purified from canine spleen.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.40 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	CA11.8H2
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab34096の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

FACS: 1/10 - 1/20.
IHC-Fr: Use at an assay dependent dilution.
IP: Use at an assay dependent dilution.

Is unsuitable for IHC-P.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性 The ad integrin is a newly recognised integrin chain that is expressed as a heterodimer with the b2 integrin (CD18). ad is expressed by a sub-population of CD8 positive T cells, macrophages and by gd T cells in the splenic red pulp.

細胞内局在 Plasma membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors