

Anti-CD134 / OX40L receptor antibody [AP-MAB0833] ab80677

製品の概要

製品名	Anti-CD134 / OX40L receptor antibody [AP-MAB0833]
製品の詳細	Rat monoclonal [AP-MAB0833] to CD134 / OX40L receptor
由来種	Rat
アプリケーション	適用あり: IHC-Fr, Flow Cyt, Cellular Activation 適用なし: Inhibition
種交差性	交差種: Mouse
免疫原	Recombinant full length protein corresponding to Mouse CD134/ OX40L receptor.
特記事項	<p>This antibody was produced from a hybridoma (mouse myeloma) fused with spleen cells from an immunized rat.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein A purified
特記事項 (精製)	The IgG fraction of culture supernatant was purified by Protein A/G affinity chromatography and lyophilized from a 0.2 µm filtered solution. Endotoxin Level: Extremely low level of LPS (< 0.002EU/µg IgG). The antibody is produced by in vitro culture.
一次抗体 備考	This antibody was produced from a hybridoma (mouse myeloma) fused with spleen cells from an immunized rat.
ポリ/モノ	モノクローナル

クローン名	AP-MAB0833
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab80677の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab18407 - Rat monoclonal IgG1, is suitable for use as an isotype control with this antibody.
Cellular Activation		Use at an assay dependent concentration. For the stimulation of T cell proliferation.

追加情報 Is unsuitable for Inhibition.

ターゲット情報

機能	Receptor for TNFSF4/OX40L/GP34.
配列類似性	Contains 4 TNFR-Cys repeats.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors