

Anti-CD30 antibody [CON6D/B5] ab23766

[3 References](#) [画像数 2](#)

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-CD30 antibody [CON6D/B5] |
| 製品の詳細 | Mouse monoclonal [CON6D/B5] to CD30 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: IHC-Fr, ELISA, IHC-P, ICC/IF |
| 種交差性 | 交差種: Human |
| 免疫原 | Hodgkin's lymphoma cell line L428 . |
| ポジティブ・コントロール | Hodgkin's lymphoma tissue or L428 cell line |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. Do Not Freeze. |
| バッファー | pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 1% BSA |
| 精製度 | Protein G purified |
| 特記事項(精製) | Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant |
| ポリ/モノ | モノクローナル |
| クローン名 | CON6D/B5 |
| ミエローマ | P3-NS1/1-Ag4-1 |
| アイソタイプ | IgG2a |

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab23766の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| IHC-Fr | | Use at an assay dependent concentration. |
| ELISA | | Use at an assay dependent concentration. |
| IHC-P | | Use at an assay dependent concentration. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. |
| ICC/IF | | 1/30. |

ターゲット情報


| | |
|-------|---|
| 機能 | Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B. |
| 配列類似性 | Contains 6 TNFR-Cys repeats. |
| 翻訳後修飾 | Phosphorylated on serine and tyrosine residues. |
| 細胞内局在 | Cytoplasm and Cell membrane. |

画像



CD30 - ab23766- Human Hodgkin's lymphoma. Formalin fixed, paraffin processed tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD30 antibody
[CON6D/B5] (ab23766)



Tissue: Hodgkin's Lymphoma
ab23766 diluted 1:30

Immunocytochemistry/ Immunofluorescence - Anti-
CD30 antibody [CON6D/B5] (ab23766)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors