

# Anti-CD5 antibody [4C7] ab851

## 2 References

### 製品の概要

製品名	Anti-CD5 antibody [4C7]
製品の詳細	Mouse monoclonal [4C7] to CD5
由来種	Mouse
アプリケーション	<b>適用あり:</b> IHC-Fr, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Fusion protein.
ポジティブ・コントロール	Mantle cell lymphoma
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.30</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: PBS, Carrier protein</p>
精製度	Van Gogh Yellow Diluent
ポリ/モノ	Tissue culture supernatant
クローン名	モノクローナル
ミエローマ	4C7
アイソタイプ	unknown
	IgG1

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab851の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		
IHC-P		

## 追加情報

This antibody may be diluted to a titer of 1:50-1:100 in an ABC method.  
We suggest an incubation period of 30-60 minutes at room temperature.  
However, depending upon the fixation conditions and the staining system employed, optimal incubation and dilutions should be determined by the user.  
Heat mediated antigen retrieval of formalin fixed paraffin embedded tissues may be required before immunostaining.

## ターゲット情報

機能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
配列類似性	Contains 3 SRCR domains.
細胞内局在	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors