

### Anti-CD5 antibody [CC17] ab34070

[1 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-CD5 antibody [CC17]
製品の詳細	Mouse monoclonal [CC17] to CD5
由来種	Mouse
特異性	ab34070 is specific for bovine CD5. This clone reacts with BoCD5.1, which is a polymorphic antigen expressed on cells of all Bos Taurus animals and a small proportion of Bos Indicus animals.
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Cow
免疫原	Tissue, cells or virus corresponding to Cow CD5. Bovine thymocytes
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	CC17
ミエローマ	NS1
アイソタイプ	IgG1

**The Abpromise guarantee** **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab34070の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

**追加情報**


FACS: 1/25 - 1/100. Use 10µl of the working dilution for 10<sup>6</sup> cells in 100ul.  
 IHC-Fr: Use at an assay dependent dilution.  
 IP: Use at an assay dependent dilution.

Not yet tested in other applications.  
 Optimal dilutions/concentrations should be determined by the end user.

**ターゲット情報**

**機能** May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.  
**配列類似性** Contains 3 SRCR domains.  
**細胞内局在** Membrane.

**画像**



Staining of bovine peripheral blood lymphocytes with Mouse anti Bovine CD5.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors