

# Anti-CD57 antibody [NK-1] - BSA and Azide free ab212405

★★★★★ [2 Abreviews](#) [画像数 1](#)

### 製品の概要

製品名	Anti-CD57 antibody [NK-1] - BSA and Azide free
製品の詳細	Mouse monoclonal [NK-1] to CD57 - BSA and Azide free
由来種	Mouse
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human <b>非交差種:</b> Rat
免疫原	Tissue, cells or virus corresponding to Human B3GAT1. Human peripheral blood mononuclear cells. Database link: <a href="#">Q9P2W7</a>
ポジティブ・コントロール	Human spleen, lymph node or tonsil tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Constituent: 100% PBS
キャリア・フリー	はい
精製度	Protein A/G purified
特記事項(精製)	ab212405 is purified from Bioreactor Concentrate by Protein A/G.
ポリ/モノ	モノクローナル
クローン名	NK-1

アイソタイプ	IgM
軽鎖の種類	kappa

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab212405の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	Use a concentration of 2 - 4 µg/ml. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. for 30 minutes at RT

## ターゲット情報

**関連性** CD57 is a carbohydrate epitope that contains a sulfoglucuronyl residue in several adhesion molecules of the nervous system. Is is created by Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 (B3GAT1) The CD57 antigen is present on subsets of T cells and natural killer cells.

## 画像



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human spleen tissue labeling B3GAT1 with ab212405 at 4µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-B3GAT1 antibody [NK-1]  
- BSA and Azide free (ab212405)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors