

Anti-CD74 antibody [LN2] ab9514

KO 評価済

★★★★★ [2 Abreviews](#) [22 References](#) [画像数 7](#)

製品の概要

製品名	Anti-CD74 antibody [LN2]
製品の詳細	Mouse monoclonal [LN2] to CD74
由来種	Mouse
アプリケーション	適用あり: Flow Cyt, WB, IHC-P, Flow Cyt (Intra)
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Human CD74. SU-DHL-4 lymphoma cells
ポジティブ・コントロール	Human tonsil normal tissue lysate - total protein: Human tonsil and Human liver tissues; Raji cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.3 Preservative: 0.1% Sodium azide Constituent: 1% BSA
精製度	Protein A/G purified
ポリ/モノ	モノクローナル
クローン名	LN2
ミエローマ	unknown
アイソタイプ	IgG1

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab9514の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/10. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
WB		Use a concentration of 5 µg/ml. Detects a band of approximately 34 kDa (predicted molecular weight: 34 kDa).
IHC-P	★★★★★ (2)	1/25 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
Flow Cyt (Intra)		Use a concentration of 1 µg/ml.

ターゲット情報

機能	Plays a critical role in MHC class II antigen processing by stabilizing peptide-free class II alpha/beta heterodimers in a complex soon after their synthesis and directing transport of the complex from the endoplasmic reticulum to the endosomal/lysosomal system where the antigen processing and binding of antigenic peptides to MHC class II takes place. Serves as cell surface receptor for the cytokine MIF.
配列類似性	Contains 1 thyroglobulin type-1 domain.
細胞内局在	Cell membrane. Endoplasmic reticulum membrane. Golgi apparatus > trans-Golgi network. Endosome. Lysosome. Transits through a number of intracellular compartments in the endocytic pathway. It can either undergo proteolysis or reach the cell membrane.

画像



Western blot - Anti-CD74 antibody [LN2] (ab9514)

All lanes : Anti-CD74 antibody [LN2] (ab9514) at 5 µg/ml

Lane 1 : Wild-type Raji cell lysate

Lane 2 : CD74 CRISPR/Cas9 edited Raji cell lysate

Lane 3 : Jurkat cell lysate

Lane 4 : HepG2 cell lysate

Lysates/proteins at 30 µg per lane.


Performed under reducing conditions.

Predicted band size: 34 kDa

Observed band size: 35 kDa


Lanes 1 - 4: Merged signal (red and green). Green - ab9514 observed at 35 kDa. Red - loading control, [ab52866](#) (Rabbit anti-alpha Tubulin antibody [EP1332Y]) observed at 55 kDa.

ab9514 was shown to react with CD74 in western blot. The band observed in CD74 CRISPR/Cas9 edited cell line [ab273378](#) (CRISPR/Cas9 edited cell lysate [ab275529](#)) below 35 kDa is likely to represent a truncated form. This has not been investigated further. Membranes were blocked in fluorescent western blot (TBS-based) blocking solution before incubation with ab9514 and [ab52866](#) (Rabbit anti-alpha Tubulin antibody [EP1332Y]) overnight at 4 °C at 5 µg/ml and a 1 in 20000 dilution respectively. Blots were incubated with Goat anti-Mouse IgG H&L (IRDye® 800CW) preabsorbed ([ab216772](#)) and Goat anti-Rabbit IgG H&L (IRDye® 680RD) preabsorbed ([ab216777](#)) secondary antibodies at 1 in 20000 dilution for 1 h at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD74 antibody [LN2] (ab9514)

Formalin-fixed, paraffin-embedded human tonsil tissue stained for CD74 using ab9514 at 1/50 dilution in immunohistochemical analysis. Antigen retrieval with citrate buffer pH 6.0.



Flow Cytometry (Intracellular) - Anti-CD74 antibody [LN2] (ab9514)

Flow cytometry overlay histogram showing wild-type Raji (green line) and CD74 knockout Raji cells (**ab273378**) stained with ab9514 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Triton X-100 for 15 min. The cells were then incubated in 1x PBS containing 10µg/ml human IgG and 10% normal goat serum to block FC receptors and non-specific protein-protein interaction followed by the antibody (ab9514) (1×10^6 in 100µl at 1 µg/ml) for 30 min at 22°C.

The secondary antibody Goat anti-mouse IgG H&L (Alexa Fluor® 488, pre-adsorbed) (**ab150117**) was used at 1/2000 for 30 min at 22°C.

Isotype control antibody was mouse IgG1κ (**ab170190**) used at the same concentration and conditions as the primary antibody (wild-type Raji cells - black line; CD74 knockout Raji cells **ab273378** - grey line). Unlabelled sample was also used as a control (this line is not shown for the purpose of simplicity).

Acquisition of >5000 events were collected using a 50 mW Blue laser (488nm) and 525/40 bandpass filter.


This antibody can also be used in Raji cells fixed with 4% formaldehyde (10 min) / permeabilized with 0.1% PBS-Triton X-100 for 15 min used under the same conditions.

Flow Cytometry - Anti-CD74 antibody [LN2] (ab9514)

Overlay histogram showing Raji cells stained with ab9514 (red line). The cells were fixed with 80% methanol (5 min) and incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions. The cells were then incubated with the

antibody (ab9514, 1/10 dilution) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) (**ab96879**) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] (**ab91353**, 2µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.

Please note that Abcam do not have any data for use of this antibody on non-fixed cells. We welcome any customer feedback.



Anti-CD74 antibody [LN2] (ab9514) at 5 µg/ml + Human tonsil normal tissue lysate - total protein (**ab29615**) at 10 µg


Secondary

Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Predicted band size: 34 kDa

Observed band size: 34 kDa

Additional bands at: 37 kDa (possible post-translational modification)



Ab9514 staining CD74 in paraffin embedded Human tonsil tissue sections by Immunohistochemistry (IHC-P).



ab9514 staining CD74 in human liver tissue by Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections). Tissue was fixed with neutral buffered formalin and a heat mediated antigen retrieval step was performed using EDTA buffer pH 9.0. Samples were then blocked with 10% serum for 10 minutes at 20°C followed by incubation with the primary antibody at a 1/75 dilution for 30 minutes at 20°C. A HRP-conjugated rat anti-mouse/rabbit polyclonal was used undiluted as the secondary antibody.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD74 antibody [LN2] (ab9514)

Image courtesy of an anonymous Abreview.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors