

### Anti-CD78 antibody [60-3G2] ab24151

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-CD78 antibody [60-3G2]   |
| 製品の詳細        | Mouse monoclonal [60-3G2] to CD78   |
| 由来種          | Mouse   |
| 特異性          | ab24151 recognises CDw78.   |
| アプリケーション     | <b>適用あり:</b> IHC-Fr, Flow Cyt   |
| 種交差性         | <b>交差種:</b> Human   |
| 免疫原          | Tissue, cells or virus corresponding to Human CD78. In vitro activated B cells.   |
| ポジティブ・コントロール | Lymphoid tissue   |
| 特記事項         | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| バッファー  | <p>pH: 7.3</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituent: PBS</p>                         |
| 精製度    | Protein G purified   |
| ポリ/モノ  | モノクローナル  |
| クローン名  | 60-3G2   |
| ミエローマ  | Sp2/0  |
| アイソタイプ | IgG1   |

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab24151の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| IHC-Fr   |           | Use at an assay dependent concentration.   |
| Flow Cyt |           | Use at an assay dependent concentration.<br><b>ab170190</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. |

## ターゲット情報

**関連性** Cdw78, also known as Ba antigen, Leu21 or LO-panB-a, is present on some immature and some mature B-cells. Its expression is greatly increased after B-cell activation in vitro. It is also found on tissue macrophages and on epithelial cells but not on T-cells, NK cells, monocytes, granulocytes, or myeloid tissues. CDw78 antigen is expressed as part of the MHC class II molecule and appears similar if not identical to the alpha beta heterodimeric MHC class II antigen.

**細胞内局在** Cell membrane

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors