

Anti-Cdc13 antibody [6F11/2] ab10873

★★★★☆ [1 Abreviews](#) [7 References](#)

製品の概要

製品名 Anti-Cdc13 antibody [6F11/2]

製品の詳細 Mouse monoclonal [6F11/2] to Cdc13

由来種 Mouse

アプリケーション **適用あり:** WB

種交差性

交差が予測される動物種: Schizosaccharomyces pombe 

免疫原 Full length native protein (purified).

特記事項

Cdc13 is a fission yeast (*S. pombe*) B-type M-phase cyclin. Cdc13 binds to Cdk1 (*cdc2*), and the resulting Cdk1-Cdc13 complex controls the G2/M transition of the cell cycle. 6F11/2 can be used for detecting Cdc13 and associated Cdk1 kinase activation.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

バッファー Preservative: 0.02% Sodium azide
Constituent: 99.98% PBS

精製度 Protein A purified

一次抗体 備考 Cdc13 is a fission yeast (*S. pombe*) B-type M-phase cyclin. Cdc13 binds to Cdk1 (*cdc2*), and the resulting Cdk1-Cdc13 complex controls the G2/M transition of the cell cycle. 6F11/2 can be used for detecting Cdc13 and associated Cdk1 kinase activation.

ポリ/モノ	モノクローナル
クローン名	6F11/2
ミエローマ	Sp2
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab10873の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 48 kDa.

ターゲット情報

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors