

Anti-CEBP Alpha/CEBPA antibody ab15048

★★★★★ [2 Abreviews](#) [8 References](#) [画像数 1](#)

製品の概要

製品名	Anti-CEBP Alpha/CEBPA antibody
製品の詳細	Rabbit polyclonal to CEBP Alpha/CEBPA
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse
免疫原	Recombinant fragment corresponding to Rat CEBP Alpha/CEBPA. Database link: NM_012524
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 99% PBS, 0.1% BSA
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab15048の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (2)	Use a concentration of 2 µg/ml. Detects a band of approximately 30 kDa (predicted molecular weight: 30,43 kDa). Detects bands of approximately 30 kDa and 43 kDa representing CEBP alpha in mouse 3T3-L1 cells.

ターゲット情報

機能	C/EBP is a DNA-binding protein that recognizes two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers.
配列類似性	Belongs to the bZIP family. C/EBP subfamily. Contains 1 bZIP domain.
細胞内局在	Nucleus.

画像



Anti-CEBP Alpha/CEBPA antibody (ab15048) + 3T3-L1 cell lysate.

Predicted band size: 30,43 kDa

Western blot of CEBP alpha using ab15048 on 3T3-L1 cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors