

Anti-CEBP Alpha/CEBPA antibody [EP709Y] - BSA and Azide free ab219726

リコンビナント RabMAb

5 References [画像数 2](#)

製品の概要

製品名	Anti-CEBP Alpha/CEBPA antibody [EP709Y] - BSA and Azide free
製品の詳細	Rabbit monoclonal [EP709Y] to CEBP Alpha/CEBPA - BSA and Azide free
由来種	Rabbit
アプリケーション	適用あり: WB 適用なし: IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	U937 cell lysate
特記事項	<p>ab219726 is the carrier-free version of ab40764.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	pH: 7.20 Constituent: PBS
キャリア・フリー	はい
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EP709Y
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab219726の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Detects a band of approximately 43 kDa (predicted molecular weight: 43 kDa).

追加情報 Is unsuitable for IP.

ターゲット情報

機能	C/EBP is a DNA-binding protein that recognizes two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers.
配列類似性	Belongs to the bZIP family. C/EBP subfamily. Contains 1 bZIP domain.
細胞内局在	Nucleus.

画像



Western blot - Anti-CEBP Alpha/CEBPA antibody [EP709Y] - BSA and Azide free (ab219726)

All lanes : Anti-CEBP Alpha/CEBPA antibody [EP709Y]

(**ab40764**) at 1/1000 dilution (Purified)

Lane 1 : U-937 (Human histiocytic lymphoma monocyte) whole cell lysate

Lane 2 : THP-1 (Human monocytic leukemia monocyte) whole cell lysate

Lane 3 : Mouse liver tissue lysate

Lane 4 : Rat liver tissue lysate

Lysates/proteins at 15 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution





Predicted band size: 43 kDa

Observed band size: 43 kDa

This data was developed using **ab40764**, the same antibody clone in a different buffer formulation.

This antibody detects high backgrounds.

Why choose a recombinant antibody?

 <p>Research with confidence Consistent and reproducible results</p>	 <p>Long-term and scalable supply Recombinant technology</p>
 <p>Success from the first experiment Confirmed specificity</p>	 <p>Ethical standards compliant Animal-free production</p>

Anti-CEBP Alpha/CEBPA antibody [EP709Y] - BSA and Azide free (ab219726)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors