

Anti-CEP170 antibody ab84545

3 References [画像数 2](#)

製品の概要

製品名	Anti-CEP170 antibody
製品の詳細	Rabbit polyclonal to CEP170
由来種	Rabbit
アプリケーション	適用あり: ICC/IF
種交差性	交差種: Human 交差が予測される動物種: Rabbit, Cow, Chimpanzee 
免疫原	Synthetic peptide within Human CEP170 aa 1500 to the C-terminus. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please contact our Scientific Support team to discuss your requirements. NP_055627.2 Database link: Q5SW79
ポジティブ・コントロール	HeLa cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.1% BSA, Tris buffered saline
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab84545の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/100 - 1/500. Permeabilization with Triton-X 100 is recommended.

ターゲット情報


関連性	CEP170, centrosomal protein 170kDa, plays a role in microtubule organization. During interphase, the encoded protein localizes to the sub-distal appendages of mature centrioles, which are microtubule-based structures thought to help organize centrosomes. During mitosis, the protein associates with spindle microtubules near the centrosomes.
細胞内局在	Centrosome, centriole. Cytoplasm, cytoskeleton. Associated with the mature mother centriole. Associated with spindle microtubules during mitosis.

画像



Immunocytochemistry/ Immunofluorescence - Anti-CEP170 antibody (ab84545)

Detection of Human CEP170 by Immunocytochemistry.
Sample: NBF-fixed asynchronous HeLa cells.
Antibody: ab84545 used at a dilution of 1/100.
Detection: Red-fluorescent Alexa Fluor 594 goat anti-rabbit IgG.



Immunocytochemistry/ Immunofluorescence - Anti-CEP170 antibody (ab84545)

Detection of Human CEP170 by Immunocytochemistry.
Sample: NBF-fixed asynchronous HeLa cells.
Antibody: ab84545 used at a dilution of 1/100.
Detection: Red-fluorescent Alexa Fluor 594 goat anti-rabbit IgG.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors