

Anti-CGA/hCG alpha antibody [INN-hFSH-132] ab11232

★★★★★ [1 Abreviews](#) [9 References](#)

製品の概要

製品名	Anti-CGA/hCG alpha antibody [INN-hFSH-132]
製品の詳細	Mouse monoclonal [INN-hFSH-132] to CGA/hCG alpha
由来種	Mouse
アプリケーション	適用あり: ELISA, IP, WB, RIA, IHC-P
種交差性	交差種: Human
免疫原	Chemical/ Small Molecule corresponding to Human CGA/hCG alpha. Human Follicle Stimulating Hormone Database link: P01215
エピトープ	ab11232 recognises the alpha 4 epitope of the alpha-subunit of Human Chorionic Gonadotrophin (hCG), hLH, hFSH, hTSH and free alpha-subunit.
特記事項	<p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein A purified
特記事項 (精製)	Purified IgG prepared by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	INN-hFSH-132

ミエローマ x63-Ag8.653

アイソタイプ IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab11232の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		
IP		
WB		
RIA		
IHC-P	★★★★★ (1)	

追加情報

ELISA: Use at an assay dependent dilution.

IHC-P: Used at a dilution of 1/1500 for 2 hrs on human placenta sections (see Abreview for further details).

IP: Use at an assay dependent dilution.

RIA: Use at an assay dependent dilution.

WB: Use at an assay dependent dilution (see Berger et al).

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

配列類似性 Belongs to the glycoprotein hormones subunit alpha family.

細胞内局在 Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors