

Anti-Choline Acetyltransferase antibody - C-terminal ab85609

★☆☆☆☆ [1 Abreviews](#) [2 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Choline Acetyltransferase antibody - C-terminal
製品の詳細	Rabbit polyclonal to Choline Acetyltransferase - C-terminal
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat
免疫原	Synthetic peptide corresponding to Rat Choline Acetyltransferase aa 600 to the C-terminus conjugated to keyhole limpet haemocyanin. (Peptide available as ab97420)
ポジティブ・コントロール	This antibody gave a positive signal in both Rat Cortex and Mouse Cholinergic nucleus Tissue Lysates.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p> <p>Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.</p>
精製度	Immunogen affinity purified

ポリ/モノ
アイソタイプ

ポリクローナル
IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab85609の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 75 kDa (predicted molecular weight: 72 kDa).

ターゲット情報

機能	Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.
関連疾患	Myasthenic syndrome, congenital, 6, presynaptic
配列類似性	Belongs to the carnitine/choline acetyltransferase family.

画像



Anti-Choline Acetyltransferase antibody - C-terminal (ab85609) at 1 µg/ml + Cortex (Rat) Tissue Lysate at 10 µg

Secondary

Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Developed using the ECL technique.


Performed under reducing conditions.

Predicted band size: 72 kDa

Observed band size: 75 kDa

Additional bands at: 18 kDa, 50 kDa. We are unsure as to the identity of these extra bands.

Exposure time: 20 minutes



Western blot - Anti-Choline Acetyltransferase antibody - C-terminal (ab85609)

Anti-Choline Acetyltransferase antibody - C-terminal (ab85609) at 1 µg/ml + Cholinergic nucleus Mouse Tissue Lysate at 10 µg

Secondary

Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) ([ab65484](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 72 kDa

Observed band size: 76 kDa

Additional bands at: 78 kDa. We are unsure as to the identity of these extra bands.

Exposure time: 2 minutes

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors