

Anti-Coronin-1C antibody ab15719

★★★★★ [1 Abreviews](#) [3 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Coronin-1C antibody
製品の詳細	Goat polyclonal to Coronin-1C
由来種	Goat
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Coronin-1C aa 450-550 (C terminal). Run BLAST with ExPASy Run BLAST with NCBI
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 0.5% BSA
精製度	Immunogen affinity purified
特記事項(精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab15719の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 0.2 - 1 µg/ml. Detects a band of approximately 55 kDa (predicted molecular weight: 55 kDa). A 1 hour primary incubation is recommended for this product.


ターゲット情報

機能	May be involved in cytokinesis, motility, and signal transduction.
組織特異性	Ubiquitous.
配列類似性	Belongs to the WD repeat coronin family. Contains 4 WD repeats.

画像



Ab15719 at 0.2µg/ml staining Coronin-1C in human lung lysate by Western blot (ECL)



Anti-Coronin-1C antibody (ab15719) at 1/1000 dilution + Mouse 3T3 fibroblast whole cell lysate (25ug)

Secondary

HRP conjugated bovine anti-goat IgG

Predicted band size: 55 kDa

Observed band size: 55 kDa

This blot shows significant non-specific background staining. Increasing the length of the blocking step and including a blocking agent in the primary antibody solution may reduce this.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors