

Anti-CSF-1-R antibody [MM0467-5V36] ab89907

1 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-CSF-1-R antibody [MM0467-5V36] |
| 製品の詳細 | Mouse monoclonal [MM0467-5V36] to CSF-1-R |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB, Neutralising |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant fragment. This information is proprietary to Abcam and/or its suppliers. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Constituent: PBS |
| 精製度 | Protein G purified |
| 特記事項(精製) | The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and 0.2 µm filtered. |
| ポリ/モノ | モノクローナル |
| クローン名 | MM0467-5V36 |
| アイソタイプ | IgG1 |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab89907の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|--------------|-----------|--|
| WB | | Use at an assay dependent concentration. |
| Neutralising | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|--|
| 機能 | Protein tyrosine-kinase transmembrane receptor for CSF1 and IL34. |
| 組織特異性 | Expressed in bone marrow and in differentiated blood mononuclear cells. |
| 配列類似性 | Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1/PDGF receptor subfamily. Contains 5 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 protein kinase domain. |
| 細胞内局在 | Membrane. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors