

Anti-CUEDC2 antibody ab79036

★★★★★ [2 Abreviews](#) [1 References](#) [画像数 1](#)

製品の概要

製品名	Anti-CUEDC2 antibody
製品の詳細	Rabbit polyclonal to CUEDC2
由来種	Rabbit
特異性	ab79036 will not cross-react with CUEDC1
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	A 19 amino acid synthetic peptide from near the carboxy terminus of human CUEDC2 (NP_076945).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Immunogen affinity purified
特記事項(精製)	ab79036 is purified by affinity chromatography via a peptide column
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab79036の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 32 kDa.

ターゲット情報

機能	Controls PGR and ESR1 protein levels through their targeting for ubiquitination and subsequent proteasomal degradation.
配列類似性	Belongs to the CUEDC2 family. Contains 1 CUE domain.
ドメイン	The CUE domain mediates interaction with PGR and ESR1.
細胞内局在	Cytoplasm. Nucleus.

画像



Lane 1 : Anti-CUEDC2 antibody (ab79036) at 1 µg/ml

Lane 2 : Anti-CUEDC2 antibody (ab79036) at 2 µg/ml

All lanes : HeLa cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 32 kDa

Observed band size: 33 kDa

Additional bands at: 50 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors