

Anti-CXCL1/2/3 (pan GRO) antibody [MM0214-11B27] ab89324

1 References [画像数 1](#)

製品の概要

製品名	Anti-CXCL1/2/3 (pan GRO) antibody [MM0214-11B27]
製品の詳細	Mouse monoclonal [MM0214-11B27] to CXCL1/GRO alpha + CXCL2 + GRO gamma
由来種	Mouse
アプリケーション	適用あり: Neutralising, WB
種交差性	交差種: Human
免疫原	Recombinant full length protein corresponding to Human GRO gamma aa 1-150. Database link: <u>P19876</u>
ポジティブ・コントロール	WB: Human placenta tissue lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein G purified
特記事項 (精製)	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	MM0214-11B27
アイソタイプ	IgG1

アプリケーション

Abpromise保証は、 次のテスト済みアプリケーションにおけるab89324の使用に適用されます

アプリケーション	Abreviews	特記事項
Neutralising		Use at an assay dependent concentration.
WB		1/500 - 1/1000. Predicted molecular weight: 11 kDa.

CXCL1/GRO alpha: Secreted. CXCL2: Secreted. GRO gamma: Secreted.

Western blot analysis of CXCL1/2/3 antibody specificity. The blot shows a single band at approximately 8 kDa in lane 1, indicated by a black arrowhead. Molecular weight markers are labeled on the left: 53, 37, 29, 20, 9, and 4 kDa.

Predicted band size: 11 kDa
Observed band size: 8 kDa

2

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors