

Anti-CYP3A5 antibody [EPR4396] ab108624

リコンビナント **RabMAb**

7 References **画像数 13**

製品の概要

製品名	Anti-CYP3A5 antibody [EPR4396]
製品の詳細	Rabbit monoclonal [EPR4396] to CYP3A5
由来種	Rabbit
アプリケーション	適用あり: WB, IP, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Human fetal liver and colon tissue lysates. IHC-P: Human liver and kidney tissues. IP: Human fetal liver tissue lysate.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Stable for 12 months at -20°C.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS, 0.05% BSA</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル

クローン名	EPR4396
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab108624の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000. Predicted molecular weight: 57 kDa. For unpurified use at 1/400.
IP		1/50. For unpurified use at 1/20.
IHC-P		1/250. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. See <u>IHC antigen retrieval protocols</u> . For unpurified use at 1/80.

ターゲット情報

機能	Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acids, and xenobiotics.
配列類似性	Belongs to the cytochrome P450 family.
細胞内局在	Endoplasmic reticulum membrane. Microsome membrane.


画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of normal human colon tissue labelling CYP3A5 with unpurified ab108624. Positive staining is shown.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody
[EPR4396] (ab108624)



Western blot - Anti-CYP3A5 antibody [EPR4396]
(ab108624)

Anti-CYP3A5 antibody [EPR4396] (ab108624) at 1/400 dilution
(unpurified) + Human fetal liver tissue lysate at 20 µg

Secondary


HRP-conjugated anti-rabbit IgG, specific to the non-reduced form of
IgG at 1/1000 dilution

Predicted band size: 57 kDa

Observed band size: 52 kDa

Blocking buffer and concentration: 5% NFDM/TBST.

Diluting buffer and concentration: 5% NFDM /TBST.




Immunoprecipitation - Anti-CYP3A5 antibody
[EPR4396] (ab108624)

ab108624 (unpurified) at 1/20 immunoprecipitating CYP3A5 in
human fetal liver tissue lysate (Lane 1). Lane 2 - PBS. For western
blotting, a HRP-conjugated anti-rabbit IgG, specific to the non-
reduced form of IgG was used as the secondary antibody (1/1500).


Blocking buffer and concentration: 5% NFDM/TBST.

Diluting buffer and concentration: 5% NFDM /TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human liver tissue labelling CYP3A5 with unpurified ab108624 at 1/80. Heat mediated antigen retrieval was performed using Tris/EDTA buffer pH 9. A prediluted HRP-polymer conjugated anti-rabbit IgG was used as the secondary antibody. Counterstained with Hematoxylin.



Western blot - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Anti-CYP3A5 antibody [EPR4396] (ab108624) at 1/1000 dilution (purified) + Human fetal liver tissue lysate at 20 µg

Secondary


HRP-conjugated anti-rabbit IgG, specific to the non-reduced form of IgG at 1/1000 dilution

Predicted band size: 57 kDa

Observed band size: 52 kDa

Blocking buffer and concentration: 5% NFDM/TBST.


Diluting buffer and concentration: 5% NFDM /TBST.



ab108624 (purified) at 1/50 immunoprecipitating CYP3A5 in human fetal liver tissue lysate (Lane 1). Lane 2 - PBS. For western blotting, a HRP-conjugated anti-rabbit IgG, specific to the non-reduced form of IgG was used as the secondary antibody (1/1500).


Blocking buffer and concentration: 5% NFDM/TBST.

Diluting buffer and concentration: 5% NFDM /TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of normal human kidney tissue labelling CYP3A5 with unpurified ab108624. Positive staining is shown.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.




All lanes : Anti-CYP3A5 antibody [EPR4396] (ab108624) at 1/1000 dilution (unpurified)

Lane 1 : Human fetal liver tissue lysate

Lane 2 : Human fetal colon tissue lysate

Lysates/proteins at 10 µg per lane.


Predicted band size: 57 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of normal human breast tissue labelling CYP3A5 with unpurified ab108624. Negative staining is shown.


Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human skeletal muscle tissue labelling CYP3A5 with unpurified ab108624. Negative staining is shown.


Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of normal human tonsil tissue labelling CYP3A5 with unpurified ab108624. Negative staining is shown.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human liver tissue labelling CYP3A5 with purified ab108624 at 1/250. Heat mediated antigen retrieval was performed using Tris/EDTA buffer pH 9. A prediluted HRP-polymer conjugated anti-rabbit IgG was used as the secondary antibody. Counterstained with Hematoxylin.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYP3A5 antibody [EPR4396] (ab108624)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-CYP3A5 antibody [EPR4396] (ab108624)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors