

Anti-Cystatin C antibody [EPR4413] ab109508

KO 評価済 リコンビナント RabMAb

★★★★★ **2 Abreviews** **24 References** 画像数 9

製品の概要

製品名	Anti-Cystatin C antibody [EPR4413]
製品の詳細	Rabbit monoclonal [EPR4413] to Cystatin C
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P 適用なし: Flow Cyt, ICC/IF or IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Human testis, Caco2, U20S, HeLa, HepG2, and U87-MG lysates, rat heart tissue lysate and mouse thymus tissue lysate. IHC-P: Human brain and kidney tissue, mouse kidney tissue and rat kidney tissue.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR4413

アイソタイプ

IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab109508の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	1/10000 - 1/50000. Detects a band of approximately 13 kDa (predicted molecular weight: 15 kDa).
IHC-P		1/250 - 1/8000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.


追加情報

Is unsuitable for Flow Cyt, ICC/IF or IP.

ターゲット情報

機能	As an inhibitor of cysteine proteinases, this protein is thought to serve an important physiological role as a local regulator of this enzyme activity.
組織特異性	Expressed in submandibular and sublingual saliva but not in parotid saliva (at protein level). Expressed in various body fluids, such as the cerebrospinal fluid and plasma. Expressed in highest levels in the epididymis, vas deferens, brain, thymus, and ovary and the lowest in the submandibular gland.
関連疾患	Amyloidosis 6 Macular degeneration, age-related, 11
配列類似性	Belongs to the cystatin family.
翻訳後修飾	The Thr-25 variant is O-glycosylated with a core 1 or possibly core 8 glycan. The signal peptide of the O-glycosylated Thr-25 variant is cleaved between Ala-20 and Val-21.
細胞内局在	Secreted.

画像



Western blot - Anti-Cystatin C antibody [EPR4413] (ab109508)

Lane 1: Wild type HAP1 whole cell lysate (20 µg)


Lane 2: CST3 knockout HAP1 whole cell lysate (20 µg)

Lane 3: Caco2 whole cell lysate (20 µg)

Lane 4: U2OS whole cell lysate (20 µg)

Lanes 1 - 4: Merged signal (red and green). Green - ab109508 observed at 15 kDa. Red - loading control, **ab9484**, observed at 37 kDa.


ab109508 was shown to specifically react with CST3 when CST3 knockout samples were used. Wild-type and CST3 knockout samples were subjected to SDS-PAGE. Ab109508 and **ab9484** (Mouse anti GAPDH loading control) were incubated overnight at 4°C at 10000 dilution and 1/10000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed **ab216773** and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed **ab216776** secondary antibodies at 1/10000 dilution for 1 hour at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cystatin C antibody [EPR4413] (ab109508)

Immunohistochemistry analysis of paraffin-embedded human kidney tissue sections labelling Cystatin C with ab109508 at 1/8000 dilution. The section was incubated with ab109508 for 30 mins at room temperature. Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**) was used as the secondary antibody. Sections were counterstained with Hematoxylin. Antigen retrieval was heat mediated with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20 mins.

Positive staining on human kidney tissue. The immunostaining was performed on a Leica Biosystems BOND® RX instrument.



Anti-Cystatin C antibody [EPR4413] (ab109508) at 1/10000 dilution
+ Mouse thymus tissue lysate at 20 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


Predicted band size: 15 kDa

Observed band size: 13 kDa

Exposure time: 3 minutes


Blocking buffer and concentration : 5% NFDM/TBST

Diluting buffer and concentration : 5% NFDM/TBST



Immunohistochemistry analysis of paraffin-embedded rat kidney tissue sections labelling Cystatin C with ab109508 at 1/8000 dilution. The section was incubated with ab109508 for 30 mins at room temperature. Rabbit specific IHC polymer detection kit HRP/DAB ([ab209101](#)) was used as the secondary antibody. Sections were counterstained with Hematoxylin. Antigen retrieval was heat mediated with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20 mins.


Positive staining on rat kidney tissue. The immunostaining was performed on a Leica Biosystems BOND® RX instrument.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cystatin C antibody [EPR4413] (ab109508)

Immunohistochemistry analysis of paraffin-embedded mouse kidney tissue sections labelling Cystatin C with ab109508 at 1/8000 dilution. The section was incubated with ab109508 for 30 mins at room temperature. Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**) was used as the secondary antibody. Sections were counterstained with Hematoxylin. Antigen retrieval was heat mediated with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20 mins.

Positive staining on mouse kidney tissue. The immunostaining was performed on a Leica Biosystems BOND® RX instrument.



Western blot - Anti-Cystatin C antibody [EPR4413] (ab109508)

Anti-Cystatin C antibody [EPR4413] (ab109508) at 1/10000 dilution + Rat heart tissue lysate at 20 µg

Secondary

Goat Anti-Rabbit IgG (HRP) with minimal cross-reactivity with human IgG at 1/2000 dilution


Predicted band size: 15 kDa

Observed band size: 13 kDa

Exposure time: 80 seconds

Blocking buffer and concentration : 5% NFDM/TBST

Diluting buffer and concentration : 5% NFDM/TBST



Western blot - Anti-Cystatin C antibody [EPR4413]
(ab109508)

All lanes : Anti-Cystatin C antibody [EPR4413] (ab109508) at
1/10000 dilution

Lane 1 : Human testis lysate

Lane 2 : HeLa cell lysate


Lane 3 : HepG2 cell lysate

Lane 4 : U87-MG cell lysate

Lysates/proteins at 10 µg per lane.

Predicted band size: 15 kDa

Observed band size: 13 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-Cystatin C antibody
[EPR4413] (ab109508)

Immunohistochemical analysis of paraffin-embedded Human brain
tissue using ab109508 at a dilution of 1/250.

Perform heat mediated antigen retrieval with citrate buffer pH 6
before commencing with IHC staining protocol.

Why choose a
recombinant antibody?



**Research with
confidence**
Consistent and
reproducible results



**Long-term and
scalable supply**
Recombinant
technology



**Success from the
first experiment**
Confirmed
specificity



**Ethical standards
compliant**
Animal-free
production

Anti-Cystatin C antibody [EPR4413] (ab109508)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors