

### Anti-Cytokeratin 5 + 6 antibody [D5/16 B4] ab17133

★★★★★ [2 Abreviews](#) [12 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-Cytokeratin 5 + 6 antibody [D5/16 B4]
製品の詳細	Mouse monoclonal [D5/16 B4] to Cytokeratin 5 + 6
由来種	Mouse
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Full length native protein (purified) corresponding to Human Cytokeratin 5 + 6.
ポジティブ・コントロール	IHC-P: Human tonsil tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.3</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituent: 1% BSA</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	D5/16 B4
アイソタイプ	IgG1
軽鎖の種類	kappa


**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab17133の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★☆ (2)	1/25 - 1/50. Pretreat tissue with citrate buffer pH 6. Incubate with primary antibody for 30 minutes at room temperature.

## ターゲット情報

関連性	Cytokeratin 5 (58 kD) is a high molecular weight, basic type of cytokeratin expressed in basal, the intermediate and the superficial cell layers of stratified epithelia as well as transitional epithelia, complex epithelia and in mesothelial cells and mesothelioma. Cytokeratin 6 (56 kD) is also a high molecular weight, basic type cytokeratin expressed by proliferating squamous epithelium often paired with Cytokeratin 16.
細胞内局在	Cytoplasmic

## 画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human tonsil tissue labeling Cytokeratin 5 + 6 with ab17133.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 5 + 6 antibody [D5/16 B4] (ab17133)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors