



Anti-DAP3 antibody [42C617] ab11928

11 References [画像数 1](#)

製品の概要

製品名	Anti-DAP3 antibody [42C617]
製品の詳細	Mouse monoclonal [42C617] to DAP3
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide:corresponding to amino acids 45-52 and 358-371 of Human DAP3.  Run BLAST with ExPASy  Run BLAST with NCBI
ポジティブ・コントロール	Jurkat whole-cell lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	42C617
アイソタイプ	IgG1

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab11928の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 45 kDa.

ターゲット情報

機能	Involved in mediating interferon-gamma-induced cell death.
組織特異性	Ubiquitous.
細胞内局在	Mitochondrion.

画像



Anti-DAP3 antibody [42C617] (ab11928) at 2 µg/ml + HepG2
(Human hepatocellular liver carcinoma cell line) Whole Cell Lysate
at 10 µg

Secondary

Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at
1/3000 dilution

Observed band size: 43 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors