

# Anti-DCV capsid polyprotein antibody ab92954

## 1 Abreviews 2 References

### 製品の概要

製品名	Anti-DCV capsid polyprotein antibody
製品の詳細	Rabbit polyclonal to DCV capsid polyprotein
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Drosophila C virus
免疫原	Synthetic peptide corresponding to C terminal residues of Drosophila C virus DCV capsid polyprotein.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 50% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
特記事項(精製)	This antibody is >90% pure
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab92954の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		1/2000 - 1/5000.

## ターゲット情報

関連性	DCV capsid polyprotein (DCVgp2) includes rhv_like or Picornavirus capsid protein domain_like and CRPV_capsid domains. Picornaviruses are non-enveloped plus-strand ssRNA animal viruses with icosahedral capsids composed of 60 copies each of 4 virus encoded proteins. Alignment includes picornaviridae, like poliovirus, hepatitis A virus, rhinovirus, foot-and-mouth disease virus and encephalomyocarditis virus. The common structure of DCV capsid polyprotein is an 8-stranded beta sandwich.
-----	---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors