

Anti-DDB2 antibody [EPR9811] ab181136

KO 評価済 リコンビナント RabMAb

9 References 画像数 5

製品の概要

製品名	Anti-DDB2 antibody [EPR9811]
製品の詳細	Rabbit monoclonal [EPR9811] to DDB2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P 適用なし: Flow Cyt or ICC/IF
種交差性	交差種: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: HeLa, SH-SY5Y and Raji cell lysates; IHC-P: Human stomach carcinoma tissue.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR9811
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab181136の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. Detects a band of approximately 45 kDa (predicted molecular weight: 48 kDa).
IHC-P		1/2500. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. See IHC antigen retrieval protocols . For unpurified use at 1/50 - 1/100.

追加情報 Is unsuitable for Flow Cyt or ICC/IF.

ターゲット情報

機能	Required for DNA repair. Binds to DDB1 to form the UV-damaged DNA-binding protein complex (the UV-DDB complex). The UV-DDB complex may recognize UV-induced DNA damage and recruit proteins of the nucleotide excision repair pathway (the NER pathway) to initiate DNA repair. The UV-DDB complex preferentially binds to cyclobutane pyrimidine dimers (CPD), 6-4 photoproducts (6-4 PP), apurinic sites and short mismatches. Also appears to function as the substrate recognition module for the DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex DDB1-CUL4-ROC1 (also known as CUL4-DDB-ROC1 and CUL4-DDB-RBX1). The DDB1-CUL4-ROC1 complex may ubiquitinate histone H2A, histone H3 and histone H4 at sites of UV-induced DNA damage. The ubiquitination of histones may facilitate their removal from the nucleosome and promote subsequent DNA repair. The DDB1-CUL4-ROC1 complex also ubiquitinates XPC, which may enhance DNA-binding by XPC and promote NER. Isoform D1 and isoform D2 inhibit UV-damaged DNA repair.
組織特異性	Ubiquitously expressed; with highest levels in corneal endothelium and lowest levels in brain. Isoform D1 is highly expressed in brain and heart. Isoform D2, isoform D3 and isoform D4 are weakly expressed.
パスウェイ	Protein modification; protein ubiquitination.
関連疾患	Defects in DDB2 are a cause of xeroderma pigmentosum complementation group E (XP-E) [MIM:278740]; also known as xeroderma pigmentosum V (XP5). XP-E is a rare human autosomal recessive disease characterized by solar sensitivity, high predisposition for developing cancers on areas exposed to sunlight and, in some cases, neurological abnormalities.
配列類似性	Belongs to the WD repeat DDB2/WDR76 family. Contains 5 WD repeats.
ドメイン	The DWD box is required for interaction with DDB1.
翻訳後修飾	Phosphorylation by ABL1 negatively regulate UV-DDB activity. Ubiquitinated by CUL4A in response to UV irradiation. Ubiquitination appears to both impair DNA-binding and promotes ubiquitin-dependent proteolysis. Degradation of DDB2 at sites of DNA damage may be a prerequisite for their recognition by XPC and subsequent repair. CUL4A-mediated degradation appears to be promoted by ABL1.

画像



Lane 1 : Anti-DDB2 antibody [EPR9811] (ab181136) at 1/2000 dilution

Lane 2 : Anti-DDB2 antibody [EPR9811] (ab181136) at 1/2000 dilution (Purified)

Lane 1 : SH-SY5Y (Human neuroblastoma epithelial cell) whole cell lysate

Lane 2 : Raji (Human Burkitt's lymphoma B lymphocyte) whole cell lysate

Lysates/proteins at 15 µg per lane.


Secondary

Lane 1 : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution (Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated)

Lane 2 : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 48 kDa

40kDa band might be an isomer of DNA damage-binding protein 2 (Isomer D3).




Anti-DDB2 antibody [EPR9811] (ab181136) at 1/1000 dilution + HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysate at 15 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) (Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated)

Predicted band size: 48 kDa

40kDa band might be an isomer of DNA damage-binding protein 2 (Isomer D3).



Western blot - Anti-DDB2 antibody [EPR9811]
(ab181136)

All lanes : Anti-DDB2 antibody [EPR9811] (ab181136) at 1/1000 dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : DDB2 knockout HeLa cell lysate

Lane 3 : Raji cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 48 kDa

Observed band size: 48 kDa

Lanes 1-3: Merged signal (red and green). Green - ab181136 observed at 48 kDa. Red - loading control, **ab8245** observed at 37 kDa.

ab181136 Anti-DDB2 antibody [EPR9811] was shown to specifically react with DDB2 in wild-type HeLa cells. Loss of signal was observed when knockout cell line **ab265115** (knockout cell lysate **ab257177**) was used. Wild-type and DDB2 knockout samples were subjected to SDS-PAGE. ab181136 and Anti-GAPDH antibody [6C5] - Loading Control (**ab8245**) were incubated overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed (**ab216776**) secondary antibodies at 1 in 10000 dilution for 1 hour at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human stomach carcinoma tissue sections labeling DDB2 with Purified ab181136 at 1/2500 dilution (0.25 µg/mL). Perform heat mediated antigen retrieval using **ab93684** (Tris/EDTA buffer, pH 9.0). ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DDB2 antibody [EPR9811] (ab181136)

Why choose a recombinant antibody?

 <p>Research with confidence Consistent and reproducible results</p>	 <p>Long-term and scalable supply Recombinant technology</p>
 <p>Success from the first experiment Confirmed specificity</p>	 <p>Ethical standards compliant Animal-free production</p>

Anti-DDB2 antibody [EPR9811] (ab181136)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors