

Anti-DeltaC antibody [zdc2] ab73336

★★★★★ **1 Abreviews** **5 References**

製品の概要

製品名	Anti-DeltaC antibody [zdc2]
製品の詳細	Mouse monoclonal [zdc2] to DeltaC
由来種	Mouse
特異性	Delta proteins activate Notch through a binding reaction that depends on their extracellular domains. During somitogenesis, Notch signalling via DeltaC keeps the oscillations of neighbouring cells in synchrony. DeltaC is also expressed in subsets of cells in the nervous system and elsewhere.
アプリケーション	適用あり: ELISA, IHC-Fr, ICC/IF, Flow Cyt, IHC - Wholemount, WB
種交差性	交差種: Zebrafish
免疫原	Extracellular region of DeltaC protein (Zebrafish).
エピトープ	DeltaC extracellular domain.
ポジティブ・コントロール	HEK293 cells transfected with DeltaC or zebrafish embryo cryosections (formulin fixed).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	zdc2

ミエローマ	Sp2/0-Ag14
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab73336の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
IHC - Wholemount	★★★★★ (1)	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

ターゲット情報

関連性 DeltaC acts as a ligand for Notch receptors and is involved in somitogenesis. It can activate Notch receptors. It is required in somite segmentation to keep the oscillations of neighboring presomitic mesoderm cells synchronized.

細胞内局在 Membrane; Single-pass type I membrane protein

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors