

Anti-Dengue Virus 1+2+3+4 antibody ab26837

★★★★★ [1 Abreviews](#) [3 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Dengue Virus 1+2+3+4 antibody
製品の詳細	Rabbit polyclonal to Dengue Virus 1+2+3+4
由来種	Rabbit
アプリケーション	適用あり: ICC/IF
種交差性	交差種: Dengue virus 1, Dengue virus 2, Dengue virus 3, Dengue virus 4
免疫原	Viral particles from C6/36 [mosquito cell line] amplified dengue type 2 virus.
ポジティブ・コントロール	dengue type III infected BHK21 cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
精製度	Whole antiserum
特記事項(精製)	Ammonium sulphate precipitated
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab26837の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (1)	

追加情報

ICC/IF: Use at a dilution of 1/50-1/100.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性

Dengue and dengue hemorrhagic fever (DHF) are caused by one of four closely related, but antigenically distinct, virus serotypes (DEN1, DEN2, DEN3, and DEN4), of the genus Flavivirus.

画像



Immunofluorescent analysis of dengue type 3 infected BHK21 cells, using 10 µg/ml ab26837. Note that not every cells are infected.

Immunocytochemistry/ Immunofluorescence - Anti-Dengue Virus 1+2+3+4 antibody (ab26837)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors