

Anti-Dexamethasone antibody ab35000

5 References

製品の概要

製品名	Anti-Dexamethasone antibody
製品の詳細	Rabbit polyclonal to Dexamethasone
由来種	Rabbit
特異性	Cross reactivity: Dexamethasone 100% Dexamethasone-21-Ac 100% Beta-Methasone 4% Dexamethasone-21-Phosphate 2% Testosterone <0.001%
アプリケーション	適用あり: ELISA
種交差性	交差種: Species independent
免疫原	Chemical/ Small Molecule: Dexamethasone conjugated to BSA.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
特記事項(精製)	This antibody was purified by ammonium sulphate fractionation.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab35000の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/1000 - 1/4000.

ターゲット情報

関連性 Dexamethasone is potent glucocorticoid anti-inflammatory agent which regulates T-cell survival, growth and differentiation. It is known by various trivial and systematic names including 9-alpha fluoro 16-alpha methylprednisolone. Dexamethasone is also used as an inducer of gene expression using a dexamethasone-inducible promoter system. In horse racing Dexamethasone is one of many steroids drugs which are prohibited but is commonly used in racing stables as a treatment for respiratory complaints.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors