

Anti-dis2 antibody ab122909

3 References [画像数 2](#)

製品の概要

製品名	Anti-dis2 antibody
製品の詳細	Rabbit polyclonal to dis2
由来種	Rabbit
特異性	ab122909 is specific to dis2 and does not cross-react with Sds21.
アプリケーション	適用あり: WB, IP
種交差性	交差種: Schizosaccharomyces pombe
免疫原	Recombinant full length Schizosaccharomyces pombe dis2.
ポジティブ・コントロール	S.pombe cell extract.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C.
バッファー	pH: 6 Preservative: 0.05% Sodium azide Constituent: 99% Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab122909の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/300 - 1/1000. Predicted molecular weight: 37 kDa.
IP		Use at an assay dependent dilution.

ターゲット情報

関連性	Essential role in cell cycle control. PP1 is perhaps required for exit from mitosis.
細胞内局在	Nucleus.

画像




All lanes : Anti-dis2 antibody (ab122909) at 1/300 dilution

Lane 1 : S.pombe cells wild type

Lane 2 : S.pombe cells dis2 deletion mutant

Predicted band size: 37 kDa



Immunoprecipitation of wild-type S. pombe extracts was performed using ab122909. Resulting supernatant (S) and immunoprecipitate (P) were immunoblotted by anti-dis2 antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors