


### Anti-DNA PKcs antibody [Y393] ab32566

KO 評価済 リコンビナント RabMAb

★★★★★ **13 Abreviews** **36 References** 画像数 12

#### 製品の概要

製品名	Anti-DNA PKcs antibody [Y393]
製品の詳細	Rabbit monoclonal [Y393] to DNA PKcs
由来種	Rabbit
特異性	This antibody is specific for DNA PKcs. It may also detect the splice isoform 2. Mouse and rat species are recommended based on WB results, we do not guarantee IHC-P for mouse and rat.
アプリケーション	<b>適用あり:</b> ChIC/CUT&RUN-seq, WB, IHC-P, ICC/IF <b>適用なし:</b> Flow Cyt or IP
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Mouse, Rat, Armenian hamster 
免疫原	Synthetic peptide within Human DNA PKcs aa 4050 to the C-terminus (C terminal). The exact sequence is proprietary. Database link: <a href="#">P78527</a>
ポジティブ・コントロール	WB: K562, MOLT4, SH-SY5Y, HeLa and Wild-type HAP1 cell lysate; Wild-type A549 cell lysate; HDLM-2 cell lysate. ICC/IF: HeLa cells. IHC-P: Human breast carcinoma tissue slides, human tonsil tissue. ChIC/CUT&RUN-Seq: U2OS cells.
特記事項	Mouse and rat samples are recommended based on WB results. This product is a recombinant monoclonal antibody, which offers several advantages including: - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information <a href="#">see here</a> . Our RabMAb <sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a> .

#### 製品の特性

製品の状態 Liquid

保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	Y393
アイソタイプ	IgG

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab32566の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ChIC/CUT&RUN-seq		Use at an assay dependent concentration. 5 µg
WB	★★★★★ (10)	1/1000 - 1/10000. Detects a band of approximately 460 kDa (predicted molecular weight: 469 kDa). <b>For unpurified, use 1/1000 - 1/2000.</b>
IHC-P		1/50. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. <b>For unpurified, use 1/5.</b>
ICC/IF	★★★★★ (2)	1/100. <b>For unpurified, use 1/10.</b>

**追加情報**      Is unsuitable for Flow Cyt or IP.

## ターゲット情報

**機能**

Serine/threonine-protein kinase that acts as a molecular sensor for DNA damage. Involved in DNA nonhomologous end joining (NHEJ) required for double-strand break (DSB) repair and V(D)J recombination. Must be bound to DNA to express its catalytic properties. Promotes processing of hairpin DNA structures in V(D)J recombination by activation of the hairpin endonuclease artemis (DCLRE1C). The assembly of the DNA-PK complex at DNA ends is also required for the NHEJ ligation step. Required to protect and align broken ends of DNA. May also act as a scaffold protein to aid the localization of DNA repair proteins to the site of damage. Found at the ends of chromosomes, suggesting a further role in the maintenance of telomeric stability and the prevention of chromosomal end fusion. Also involved in modulation of transcription. Recognizes the substrate consensus sequence [ST]-Q. Phosphorylates 'Ser-139' of histone variant H2AX/H2AFX, thereby regulating DNA damage response mechanism. Phosphorylates DCLRE1C, c-Abl/ABL1, histone H1, HSPCA, c-jun/JUN, p53/TP53, PARP1, POU2F1, DHX9, SRF, XRCC1, XRCC1, XRCC4, XRCC5, XRCC6, WRN, MYC and RFA2. Can phosphorylate C1D not only in the presence of linear DNA but also in the presence of supercoiled DNA. Ability to phosphorylate p53/TP53 in the presence of supercoiled DNA is dependent on

C1D.

## 配列類似性

Belongs to the PI3/PI4-kinase family.  
Contains 1 FAT domain.  
Contains 1 FATC domain.  
Contains 2 HEAT repeats.  
Contains 1 PI3K/PI4K domain.  
Contains 3 TPR repeats.


## 翻訳後修飾

Phosphorylated upon DNA damage, probably by ATM or ATR. Autophosphorylated on Thr-2609, Thr-2638 and Thr-2647. Thr-2609 is a DNA damage-inducible phosphorylation site (inducible with ionizing radiation, IR). Autophosphorylation induces a conformational change that leads to remodeling of the DNA-PK complex, requisite for efficient end processing and DNA repair.  
S-nitrosylated by GAPDH.

## 細胞内局在

Nucleus.


## 画像



ChIC/CUT&RUN was performed using a pAG-MNase at a final concentration of 700 ng/mL,  $2 \times 10^5$  U2OS cells and 5  $\mu$ g of ab32566 [Y393]. The resulting DNA was sequenced on the Illumina NovaSeq 6000 to a depth of 10 million reads. The negative IgG control **ab172730** is also shown.

Additional screenshots of mapped reads can be downloaded [here](#).

The University of Geneva owns patents relevant to ChIC (Chromatin Immuno-Cleavage) methods.



**All lanes** : Anti-DNA PKcs antibody [Y393] (ab32566) at 1/1000 dilution

**Lane 1** : Wild-type A549 cell lysate

**Lane 2** : PRKDC knockout A549 cell lysate

**Lane 3** : K562 cell lysate

**Lane 4** : HDLM-2 cell lysate


Lysates/proteins at 20  $\mu$ g per lane.

Performed under reducing conditions.

**Predicted band size:** 469 kDa

**Observed band size:** 450 kDa

False colour image of Western blot: Anti-DNA PKcs antibody [Y393] staining at 1/1000 dilution, shown in green; Mouse anti-CANX [CANX/1543] ([ab238078](#)) loading control staining at 1/20000 dilution, shown in red. In Western blot, ab32566 was shown to bind specifically to DNA PKcs. A band was observed at 450 kDa in wild-type A549 cell lysates with no signal observed at this size in PRKDC knockout cell line. To generate this image, wild-type and PRKDC knockout A549 cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 3 % milk in TBS-0.1 % Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4 °C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were Goat anti-Rabbit IgG H&L 800CW and Goat anti-Mouse IgG H&L 680RD at 1/20000 dilution.



**Lane 1:** Wild-type HAP1 cell lysate (20 µg)


**Lane 2:** DNA PKcs knockout HAP1 cell lysate (20 µg)

**Lane 3:** K562 cell lysate (20 µg)

**Lane 4:** SH-SY5Y cell lysate (20 µg)


**Lanes 1 - 4:** Merged signal (red and green). Green - ab32566 observed at 470 kDa. Red - loading control, [ab7291](#), observed at 52 kDa.

ab32566 was shown to specifically react with DNA PKCs in wild-type HAP1 cells. No band was observed when DNA PKcs knockout samples were examined. Wild-type and DNA PKCs knockout samples were subjected to SDS-PAGE. ab32566 and [ab7291](#) (loading control to alpha tubulin) were diluted 1/1000 and 1/10,000 and incubated overnight at 4°C. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1/10,000 dilution for 1 hour at room temperature before imaging.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DNA PKcs antibody [Y393] (ab32566)

Immunohistochemical staining of paraffin embedded human tonsil with purified ab32566 at a working dilution of 1 in 50. The secondary antibody used is a HRP polymer for rabbit IgG. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.



Immunocytochemistry/ Immunofluorescence - Anti-DNA PKcs antibody [Y393] (ab32566)

Immunofluorescent staining of HeLa cells (fixed with 4% PFA) with purified ab32566 at a dilution of 1/100. An Alexa Fluor® 555 goat anti-rabbit antibody was used as the secondary at a dilution of 1/200. The panel on the right shows the DAPI counter-staining.



Western blot - Anti-DNA PKcs antibody [Y393] (ab32566)

**All lanes :** Anti-DNA PKcs antibody [Y393] (ab32566) at 1/7300 dilution (purified)

**Lane 1 :** MOLT4 cell lysate

**Lane 2 :** HeLa cell lysate

Lysates/proteins at 20 µg per lane.

**Secondary**


**All lanes :** HRP goat anti-rabbit (H+L) at 1/1000 dilution

**Predicted band size:** 469 kDa

**Observed band size:** 460 kDa


Blocking buffer: 5% NFDm/TBST

Dilution buffer: 5% NFDm/TBST




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DNA PKcs antibody [Y393] (ab32566)

Immunohistochemical staining of paraffin embedded human tonsil with unpurified ab32566 at a working dilution of 1 in 5. The secondary antibody used is a HRP polymer for rabbit IgG. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.



Immunocytochemistry/ Immunofluorescence - Anti-DNA PKcs antibody [Y393] (ab32566)

Immunofluorescent staining of HeLa cells (fixed with 4% PFA) with unpurified ab32566 at a dilution of 1/10. An Alexa Fluor<sup>®</sup> 555 goat anti-rabbit antibody was used as the secondary at a dilution of 1/200. The panel on the right shows the DAPI counter-staining.



Western blot - Anti-DNA PKcs antibody [Y393] (ab32566)

**All lanes** : Anti-DNA PKcs antibody [Y393] (ab32566) at 1/1000 dilution (unpurified)

**Lane 1** : MOLT4 cell lysate

**Lane 2** : HeLa cell lysate

Lysates/proteins at 20 µg per lane.

**Secondary**

**All lanes** : HRP goat anti-rabbit (H+L) at 1/1000 dilution


**Predicted band size:** 469 kDa

**Observed band size:** 460 kDa

Blocking buffer: 5% NFDM/TBST


Dilution buffer: 5% NFDM/TBST





Unpurified ab32566, at a 1/50 dilution, staining DNA PKCs in paraffin embedded human breast carcinoma tissue by immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DNA PKcs antibody [Y393] (ab32566)



Anti-DNA PKcs antibody [Y393] (ab32566) at 1/1000 dilution (unpurified) + K562 cell lysate.

**Predicted band size:** 469 kDa  
**Observed band size:** 460 kDa

Western blot - Anti-DNA PKcs antibody [Y393] (ab32566)

Why choose a recombinant antibody?

 <b>Research with confidence</b> Consistent and reproducible results	 <b>Long-term and scalable supply</b> Recombinant technology
 <b>Success from the first experiment</b> Confirmed specificity	 <b>Ethical standards compliant</b> Animal-free production

Anti-DNA PKcs antibody [Y393] (ab32566)



## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors