

Anti-DPL antibody ab23701

★☆☆☆☆ [2 Abreviews](#) [画像数 1](#)

製品の概要

| | |
|----------|--|
| 製品名 | Anti-DPL antibody |
| 製品の詳細 | Rabbit polyclonal to DPL |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Mouse |
| 免疫原 | Synthetic peptide corresponding to Human DPL aa 112-120. Sequence: TQAANQGEF |

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | pH: 7.40 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 50% Glycerol, 0.1% BSA |
| 精製度 | Immunogen affinity purified |
| ポリモノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab23701の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | ★☆☆☆☆ (2) | Use a concentration of 2 µg/ml. Predicted molecular weight: 20 kDa. |

ターゲット情報

| | |
|-------|------------------------------|
| 組織特異性 | Expressed in testis. |
| 配列類似性 | Belongs to the prion family. |
| 細胞内局在 | Cell membrane. |

画像



All lanes : Anti-DPL antibody (ab23701) at 2 µg/ml

Lane 1 : 27ug mouse spleen 100k x g supernatant.

Lane 2 : 54ug mouse spleen 100k x g supernatant.

Performed under reducing conditions.

Predicted band size: 20 kDa

Observed band size: 19 kDa

Additional bands at: 25 kDa (possible glycosylated form)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors