

DyLight® 488 Anti-S tag antibody ab117509

1 References

製品の概要

製品名	DyLight® 488 Anti-S tag antibody
製品の詳細	DyLight® 488 Rabbit polyclonal to S tag
由来種	Rabbit
標識	DyLight® 488. Ex: 493nm, Em: 518nm
アプリケーション	適用あり: ICC/IF, Flow Cyt
種交差性	交差種: Species independent
免疫原	Synthetic peptide: KETAAAKFERQHMS

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.20 Constituent: 99% Tris buffered saline
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション
